

# Wheelchair Safety is Key!



*Use proper  
securement at  
all times.*

## **5** Steps:

### **Step 1**

#### **Securement Locations**

- Locate mobility device in center of the 4 floor attachments. Set wheel locks/turn off power.
- Position rear belts parallel with rear wheels.
- Front belts to flair out slightly for lateral stability.
- Belts should not bend around any part of wheelchair/scooter.

### **Step 2**

#### **Front Securement**

- Apply these first if they are of the cam locking type.
- Attach belt to frame of wheelchair at junction.

### Step 3

#### Rear Securement

- Apply to chair at a solid junction of frame, if possible about 2" below seat level.
- Straps should be about 30-45 degrees.
- Tighten securely.

### Step 4

#### Lap Belt

- Locate lap belt low at pelvic area, not over abdomen.
- Position 45-75 degrees from horizontal, snug but not tight.

### Step 5

#### Shoulder Belt

- Shoulder belt should be located attached to a lap belt.
- Locate belt over shoulder.
- If not spring tensioned, leave some slack-fist test.

